



# The Aer X Difference

It's what's on  
the inside  
that counts.

The Aer X portable oxygen concentrator by 3B™ Medical brings modern design and a smaller footprint for patients on the go.

## ACTIVE X

Turbo for your Respiratory Engine



## Introducing Active X

At 3B™ Medical we know that being on oxygen isn't limiting your desire to stay active. Life doesn't stop because of a diagnosis. For you we have created Active X technology. Active X keeps up with your oxygen requirements based on your breathing rates. During activities that increase your heart rate, such as brisk walking or climbing stairs, your body requires more oxygen and you breathe faster. Aer X uses Active X technology to detect your breaths per minute and will provide extra oxygen to fuel your Respiratory Engine.

Life doesn't stop because of a diagnosis.

It is sort of like having an on-demand turbo function for your respiratory engine.

## Magnetic Port

The Aer X's Magnetic Port is designed for safety. In the event your oxygen cannula is caught on an object, the unique Magnetic Port quickly releases from the device.





## Easy Change Sieve Cartridge

More than just good looks. With its easy slide and snap function, maintenance is a breeze!

## Oxygen Bolus

The Aer X's Constant Pulse Volume delivers up to 1,200 ml per minute, almost 20% more oxygen than competing units. This higher capacity means that patients are more comfortable ambulating with the Aer X's increased oxygen delivery.

## Proudly Made in the USA

The Aer X, by 3B™ Medical is designed and manufactured in the United States. Help your patients explore the freedom and mobility that the Aer X allows, at a value/cost point that will have you viewing POCs differently.



### Aer X Starter Pack includes:

- POC with Power Adapter AC Wall Plug and Power Cord
- Custom Carry Bag with Adjustable Strap
- Car Charger
- Rechargeable Battery
- 3 ft. Ultra Soft Cannula

### Device Specs:

Model Number:	AX1000
Height:	8.25
Depth:	2.75
Length:	7.12
Weight:	4.2 lb
Hours of Use:	4 hr
Pulse Settings:	1-5
Noise Level:	39DbA at 21lpm
Bolus:	maximum O2 per min. is 1200 ml/ min



Sales and Information contact [info@3Bproducts.com](mailto:info@3Bproducts.com) | [myAerX.com](http://myAerX.com)